

Product brochure

08.08.2026

KERN 317-009-700

POM-Box for individual weights (Fracture gram, 1 mg - 500 mg (einzel))

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



APPLICATION EXAMPLE



FEATURES

- Clean room box made of polyoxymethylene
- Especially suitable for use in the pharmaceutical industry
- Box can be disinfected
- Heat-resistant material
- Unlined

PICTOGRAMS

Standard



SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day