

KERN 317-02

Individual weight, OIML Class E2, 2 g, Knob, stainless steel polished



Article number for guaranteed plus-deviation	921-332
--	---------

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	6×6×10 mm
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Optional cases

Case - Plastic article number	317-020-400
Case - Aluminium article number	317-020-600
Case - Wood article number	317-020-100

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	317

Measuring System

Conformity assessment possible	✓
Test weight - volume (cm ³)	0,25
Test weight - volume uncertainty (cm ³)	0,0031
Test weight - Density (kg/m ³)	7950
Test weight - Density uncertainty (kg/m ³)	140
Test weight - Density [Min] (OIML)	6000
Test weight - Density [Max] (OIML)	12000
Test weight - nominal value	2 g
OIML-compliant	✓
Weight OIML class	E2
Permissible tolerance according to OIML (+/-)	0,04 mg

Construction

Shape weight	Knob
Dimensions (Ø×H)	6×10 mm
Material	stainless steel polished

Services (optional)

Article number for conformity assessment	952-332
Article number for DAkkS calibration	962-332