

Product brochure

08.08.2026

KERN 317-07

Individual weight, OIML Class E2, 100 g, Knob, stainless steel polished

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
315-243	Tweezers 315-243 (stainless steel 105 mm)
315-250	Tweezers 315-250 (stainless steel 150 mm)
317-070-100	Wood weight case 317-070-100 (for nominal values 100 g for classes E1+E2+F1 for design Cylindrical/compact)
317-070-400	Plastic weight case 317-070-400 (for nominal values 100 g for classes E1, E2 for design Cylindrical/compact)
317-070-600	Aluminium weight case 317-070-600 (for nominal values 100 g for classes E1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
921-337	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-337
962-337	DAkKS calibration certificate 962-337
962-337R	DAkKS Recalibration 962-337R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days