

Product brochure

08.08.2026

KERN 317-140-700

POM-Box for individual weights (Cylindrical/compact, 10 kg)

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



FEATURES

- Clean room box made of polyoxymethylene
- Especially suitable for use in the pharmaceutical industry
- Box can be disinfected
- Heat-resistant material
- Unlined

PICTOGRAMS

Standard



SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day