


# KERN 318-01

**KERN**

Milligram weights, E2, 1 mg, Platelet, stainless steel polished



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	318

## Measuring System

Conformity assessment possible	✓
Test weight - volume (cm <sup>3</sup> )	0,0004
Test weight - volume uncertainty (cm <sup>3</sup> )	0,0001
Test weight - Density (kg/m <sup>3</sup> )	7960
Test weight - Density uncertainty (kg/m <sup>3</sup> )	170
Test weight - nominal value	1 mg
OIML-compliant	✓
Weight OIML class	E2
Permissible tolerance according to OIML (+/-)	0,006 mg

## Construction

Shape weight	Platelet
Dimensions weight (W×D×H)	3,5×3,5×3 mm
Material	stainless steel polished

## Services (optional)

Article number for conformity assessment	952-351
Article number for DAkkS calibration	962-351
Article number for guaranteed plus-deviation	921-351

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	3,5×3,5×3 mm
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

## Optional cases

Case - Plastic article number	347-009-400
Case - Aluminium article number	317-009-600
Case - Wood article number	338-090-200

## Pictograms

### STANDARD

