

Product brochure

08.08.2026

KERN 318-42

Set of milligram weights, E2, 1 mg - 500 mg, Wire, stainless steel polished, in plastic case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
315-243	Tweezers 315-243 (stainless steel 105 mm)
315-247	Tweezers 315-247 (stainless steel, plastic 136 mm)
317-009-600	Aluminium weight case 317-009-600 (for nominal values 1 mg - 500 mg (einzel) for classes E1 - M3 for design Fracture gram)
318-270	Brush 318-270
338-090-200	Wood weight case 338-090-200 (for nominal values 1 mg - 500 mg (einzel) for classes E1 - M2 for design Cylindrical/compact)
347-009-400	Plastic weight case 347-009-400 (for nominal values 1 mg - 500 mg (einzel) for classes E1 - M2 for design Cylindrical/compact)

SERVICES

Model	Description
962-250R	DAkkS Recalibration 962-250R
962-350	DAkkS calibration certificate 962-350

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days