

Individual weight, OIML Class F1, 200 g, Eco Shape, stainless steel polished



Date of creation: 15.08.2026

Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	326

Measuring System

Test weight - nominal value	200 g
Test weight - volume uncertainty (cm ³)	0,3125
Weight OIML class	F1
Test weight - volume (cm ³)	25
Permissible tolerance according to OIML (+/-)	1 mg
OIML-compliant	✓
Test weight - Density (kg/m ³)	7960
Test weight - Density uncertainty (kg/m ³)	100

Construction

Shape weight	Eco Shape
Dimensions (Ø×H)	28×48,5 mm
Material	stainless steel polished
Adjustment method	Adjusting cavity - Seal

Services (optional)

Article number for DAkkS calibration	962-438
--------------------------------------	---------

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	28×28×47 mm
Shipping method	Parcel service
Net weight approx.	0,20 kg
Gross weight approx.	0,20 kg
Shipping weight	0,2 kg
Shipping volume weight	0,1 kg

Optional cases

Case - Plastic article number	347-080-400
Case - Aluminium article number	317-080-600
Case - Wood article number	317-080-100

Product information

GTIN/EAN number	4045761247291
Trade number	84239090
Warranty Time	24 mon
REACH classification	no prohibited substances

Pictograms

STANDARD



OPTION

