

KERN 327-05

KERN

Individual weight, OIML Class F1, 20 g, Knob, stainless steel polished



Date of creation: 15.08.2026

Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	327

Measuring System

Test weight - nominal value	20 g
Test weight - volume uncertainty (cm ³)	0,0313
Weight OIML class	F1
Test weight - volume (cm ³)	2,5
Permissible tolerance according to OIML (+/-)	0,25 mg
OIML-compliant	✓
Test weight - Density (kg/m ³)	7960
Test weight - Density uncertainty (kg/m ³)	100

Construction

Shape weight	Knob
Dimensions (Ø×H)	13×23,5 mm
Material	stainless steel polished
Adjustment method	Polishing

Services (optional)

Article number for DAkkS calibration	962-435
--------------------------------------	---------

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	13×13×23 mm
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg
Shipping volume weight	0,1 kg

Optional cases

Case - Plastic article number	347-050-400
Case - Aluminium article number	317-050-600
Case - Wood article number	317-050-100

Product information

GTIN/EAN number	4045761247116
Trade number	84239090
Warranty Time	24 mon
REACH classification	no prohibited substances

Pictograms

STANDARD



OPTION

