


Product brochure

08.08.2026

KERN 328-01

Milligram weights, F1, 1 mg, Platelet, stainless steel polished

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
315-243	Tweezers 315-243 (stainless steel 105 mm)
315-247	Tweezers 315-247 (stainless steel, plastic 136 mm)
317-009-600	Aluminium weight case 317-009-600 (for nominal values 1 mg - 500 mg (einzel) for classes E1 - M3 for design Fracture gram)
318-270	Brush 318-270
338-090-200	Wood weight case 338-090-200 (for nominal values 1 mg - 500 mg (einzel) for classes E1 - M2 for design Cylindrical/compact)
347-009-400	Plastic weight case 347-009-400 (for nominal values 1 mg - 500 mg (einzel) for classes E1 - M2 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
921-451	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-451
962-451	DAkkS calibration certificate 962-451
962-451R	DAkkS Recalibration 962-451R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days