


# KERN 328-09

# KERN

Milligram weights, F1, 500 mg, Platelet, stainless steel polished (OIML)



## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	12×12×4 mm
Net weight	0,0005 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

## Optional cases

Case - Plastic article number	347-009-400
Case - Aluminium article number	317-009-600
Case - Wood article number	338-090-200

## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	328

## Measuring System

Conformity assessment possible	✓
Test weight - volume (cm <sup>3</sup> )	0,0581
Test weight - volume uncertainty (cm <sup>3</sup> )	0,0011
Test weight - Density (kg/m <sup>3</sup> )	7960
Test weight - Density uncertainty (kg/m <sup>3</sup> )	170
Test weight - nominal value	500 mg
OIML-compliant	✓
Weight OIML class	F1
Permissible tolerance according to OIML (+/-)	0,08 mg

## Construction

Shape weight	Platelet
Dimensions weight (W×D×H)	12×12×4 mm
Material	stainless steel polished (OIML)

## Services (optional)

Article number for conformity assessment	952-459
Article number for DAkkS calibration	962-459
Article number for guaranteed plus-deviation	921-459