

Product brochure

08.08.2026

KERN 328-22

Set of milligram weights, F1, 1 mg - 500 mg, Platelet, stainless steel polished, in plastic case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
315-243	Tweezers 315-243 (stainless steel 105 mm)
315-247	Tweezers 315-247 (stainless steel, plastic 136 mm)
317-009-600	Aluminium weight case 317-009-600 (for nominal values 1 mg - 500 mg (einzel) for classes E1 - M3 for design Fracture gram)
318-270	Brush 318-270
328-01	Test weight 328-01 (OIML Class F1, Nominal value 1 mg, Platelet, stainless steel polished)
328-02	Test weight 328-02 (OIML Class F1, Nominal value 2 mg, Platelet, stainless steel polished)
328-03	Test weight 328-03 (OIML Class F1, Nominal value 5 mg, Platelet, stainless steel polished)
328-04	Test weight 328-04 (OIML Class F1, Nominal value 10 mg, Platelet, stainless steel polished)
328-05	Test weight 328-05 (OIML Class F1, Nominal value 20 mg, Platelet, stainless steel polished)
328-06	Test weight 328-06 (OIML Class F1, Nominal value 50 mg, Platelet, stainless steel polished)
328-07	Test weight 328-07 (OIML Class F1, Nominal value 100 mg, Platelet, stainless steel polished)
328-08	Test weight 328-08 (OIML Class F1, Nominal value 200 mg, Platelet, stainless steel polished)
328-09	Test weight 328-09 (OIML Class F1, Nominal value 500 mg, Platelet, stainless steel polished)

Model	Description
328-222-420	Weight case 328-222-420 (for nominal values 1 mg - 5 g for classes E1 - M2 for design Cylindrical/compact)
338-090-200	Wood weight case 338-090-200 (for nominal values 1 mg - 500 mg (einzel) for classes E1 - M2 for design Cylindrical/compact)
347-009-400	Plastic weight case 347-009-400 (for nominal values 1 mg - 500 mg (einzel) for classes E1 - M2 for design Cylindrical/compact)

SERVICES

Model	Description
920-701	Marking (Last digits of serial number) KERN 920-701
962-450	DAkKS calibration certificate 962-450
962-450R	DAkKS Recalibration 962-450R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days