

Product brochure

08.08.2026

# KERN 328-37

Milligram weights, F1, 100 mg, Wire, stainless steel polished

ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

+49 7433 9933-0  
info@kern-sohn.com  
www.kern-sohn.com

## PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

## PICTOGRAMS

Standard

Option



## ACCESSORIES

| Model       | Description                                                                                                                        |
|-------------|------------------------------------------------------------------------------------------------------------------------------------|
| 315-242     | Tweezers 315-242 (polyamide   100 mm)                                                                                              |
| 315-243     | Tweezers 315-243 (stainless steel   105 mm)                                                                                        |
| 315-247     | Tweezers 315-247 (stainless steel, plastic   136 mm)                                                                               |
| 317-009-600 | Aluminium weight case 317-009-600 (for nominal values 1 mg - 500 mg (einzel)   for classes E1 - M3   for design Fracture gram)     |
| 318-270     | Brush 318-270                                                                                                                      |
| 338-090-200 | Wood weight case 338-090-200 (for nominal values 1 mg - 500 mg (einzel)   for classes E1 - M2   for design Cylindrical/compact)    |
| 347-009-400 | Plastic weight case 347-009-400 (for nominal values 1 mg - 500 mg (einzel)   for classes E1 - M2   for design Cylindrical/compact) |

## SERVICES

| Model    | Description                           |
|----------|---------------------------------------|
| 962-457  | DAkKS calibration certificate 962-457 |
| 962-457R | DAkKS Recalibration 962-457R          |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days