

KERN 346-79



Individual weight, OIML Class M1, 50 kg, Block Eco, cast iron lacquered



Date of creation: 26.09.2026

Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	346

Measuring System

Test weight - nominal value	50 kg
Test weight - volume uncertainty (cm ³)	595,12
Weight OIML class	M1
Test weight - volume (cm ³)	7042,2535
Permissible tolerance according to OIML (+/-)	2,5 g
OIML-compliant	✓
Test weight - Density (kg/m ³)	7100
Test weight - Density uncertainty (kg/m ³)	600

Construction

Shape weight	Block Eco
Dimensions weight (W×D×H)	310×157×195 mm
Material	cast iron lacquered
Adjustment method	Adjusting cavity - Seal

Services (optional)

Article number for DAkkS calibration	962-646
--------------------------------------	---------

Packing & Shipping

Delivery time	2 d
Dimensions packaging (W×D×H)	330×180×230 mm
Shipping method	Freight
Net weight approx.	50 kg
Gross weight approx.	50 kg
Shipping weight	50 kg
Shipping volume weight	2,7 kg

Product information

Warranty Time	24 mon
REACH classification	no prohibited substances

Pictograms

STANDARD



OPTION

