


KERN 346-88

Individual weight, OIML Class M1, 20 kg, Block, cast iron lacquered



Dimensions packaging (W×D×H)	260×127×150 mm
Shipping method	Parcel service
Net weight approx.	20 kg
Gross weight approx.	22 kg
Shipping weight	20,15 kg

Optional cases

Case - Aluminium article number	346-080-600
---------------------------------	-------------

Product information

GTIN/EAN number	4045761361805
Trade number	84239090
REACH classification	no prohibited substances

Pictograms

STANDARD



OPTION



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	346

Measuring System

Test weight - nominal value	20 kg
Test weight - volume uncertainty (cm ³)	238,048
Weight OIML class	M1
Test weight - volume (cm ³)	2816,9014
Permissible tolerance according to OIML (+/-)	1 g
OIML-compliant	✓
Test weight - Density [Min] (OIML)	3000
Test weight - Density (kg/m ³)	7100
Test weight - Density uncertainty (kg/m ³)	600

Construction

Shape weight	Block
Dimensions weight (W×D×H)	230×113×140 mm
Material	cast iron lacquered
Adjustment method	Adjusting cavity - Seal

Services (optional)

Article number for DAkkS calibration	962-645
--------------------------------------	---------

Packing & Shipping

Delivery time	1 d
---------------	-----