

Product brochure

08.08.2026

KERN 347-016

Individual weight, OIML Class M1, 1 g, Hook weight, stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PICTOGRAMS

Standard

Option





ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
347-030-400	Plastic weight case 347-030-400 (for nominal values 1 g - 5 g for classes F1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-631	DAkkS calibration certificate 962-631
962-631R	DAkkS Recalibration 962-631R

ENGLISH

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days