

Product brochure

08.08.2026

KERN 347-046

Individual weight, OIML Class M1, 10 g, Hook weight, stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
347-050-400	Plastic weight case 347-050-400 (for nominal values 10 g - 20 g for classes F1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-634	DAkkS calibration certificate 962-634
962-634R	DAkkS Recalibration 962-634R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days