

# KERN 347-047

**KERN**

Individual weight, OIML Class M1, 10 g, Slot weight (Type B), stainless steelfine turned



Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

## Measuring System

Test weight - volume (cm <sup>3</sup> )	1,1905
Test weight - volume uncertainty (cm <sup>3</sup> )	0,0241
Test weight - Density (kg/m <sup>3</sup> )	8000
Test weight - Density uncertainty (kg/m <sup>3</sup> )	100
Test weight - Density [Min] (OIML)	2600
Test weight - nominal value	10 g
OIML-compliant	✓
Weight OIML class	M1
Permissible tolerance according to OIML (+/-)	2 mg

## Construction

Shape weight	Slot weight (Type B)
Dimensions (Ø×H)	28×5,8 mm
Material	stainless steelfine turned

## Services (optional)

Article number for DAkkS calibration	962-634
--------------------------------------	---------

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	28×28×6 mm
Net weight	0,01 kg
Shipping method	Parcel service