

Product brochure

08.08.2026

KERN 347-055

Individual weight, OIML Class M1, 20 g, Slot weight (Type A), stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- Compatible with all slot weights/supporting bars of type A

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
347-075-100	Supporting bar 347-075-100, OIML Class M1, Nominal value 100 g, stainless steel fine turned
347-080-400	Plastic weight case 347-080-400 (for nominal values 200 g for classes F1 - M3 for design Cylindrical/compact)
347-095-100	Supporting bar 347-095-100, OIML Class M1, Nominal value 500 g, stainless steel fine turned
347-115-100	Supporting bar 347-115-100, OIML Class M1, Nominal value 1 kg, stainless steel fine turned

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-635	DAkks calibration certificate 962-635
962-635R	DAkks Recalibration 962-635R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days