

Product brochure

08.08.2026

# KERN 347-056

Individual weight, OIML Class M1, 20 g, Hook weight, stainless steel fine turned

ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

☎ +49 7433 9933-0  
✉ [info@kern-sohn.com](mailto:info@kern-sohn.com)  
🌐 [www.kern-sohn.com](http://www.kern-sohn.com)

PICTOGRAMS

Standard

Option





ACCESSORIES

| Model       | Description                                                                                                             |
|-------------|-------------------------------------------------------------------------------------------------------------------------|
| 315-242     | Tweezers 315-242 (polyamide   100 mm)                                                                                   |
| 347-050-400 | Plastic weight case 347-050-400 (for nominal values 10 g - 20 g   for classes F1 - M3   for design Cylindrical/compact) |

SERVICES

| Model    | Description                                                         |
|----------|---------------------------------------------------------------------|
| 920-350  | Marking by laser engraving (individual marking) KERN 920-350        |
| 920-351  | Marking of test weights (last digits of serial number) KERN 920-351 |
| 962-635  | DAkkS calibration certificate 962-635                               |
| 962-635R | DAkkS Recalibration 962-635R                                        |

ENGLISH

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days