

KERN 347-065

KERN

Individual weight, OIML Class M1, 50 g, Slot weight (Type A), stainless steelfine turned



Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Optional cases

Case - Plastic article number	347-080-400
-------------------------------	-------------

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

Measuring System

Test weight - volume (cm ³)	5,9524
Test weight - volume uncertainty (cm ³)	0,1205
Test weight - Density (kg/m ³)	8000
Test weight - Density uncertainty (kg/m ³)	100
Test weight - Density [Min] (OIML)	2600
Test weight - nominal value	50 g
OIML-compliant	✓
Weight OIML class	M1
Permissible tolerance according to OIML (+/-)	3 mg

Construction

Shape weight	Slot weight (Type A)
Dimensions (Ø×H)	28×13 mm
Material	stainless steelfine turned

Services (optional)

Article number for DAkkS calibration	962-636
--------------------------------------	---------

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	28×28×13 mm
Shipping method	Parcel service