

Product brochure

08.08.2026

KERN 347-067

Individual weight, OIML Class M1, 50 g, Slot weight (Type B), stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- Compatible with all slot weights/supporting bars of type B

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
347-117-100	Supporting bar 347-117-100, OIML Class M1, Nominal value 1 kg, stainless steel fine turned
347-137-100	Supporting bar 347-137-100, OIML Class M1, Nominal value 5 kg, stainless steel fine turned
347-147-100	Supporting bar 347-147-100, OIML Class M1, Nominal value 10 kg, stainless steel fine turned

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-636	DAkkS calibration certificate 962-636
962-636R	DAkkS Recalibration 962-636R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days