

Product brochure

08.08.2026

KERN 347-085-N00

Individual weight (Newton), 2 N, Slot weight (Type A), stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- Compatible with all slot weights/supporting bars of type A

PICTOGRAMS

Standard



Option



ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
347-075-100-N00	Supporting Bar (Newton) 347-075-100-N00, Nominal value 1 N, stainless steel fine turned
347-090-400	Plastic weight case 347-090-400 (for nominal values 500 g for classes F1 - M3 for design Cylindrical/compact)
347-095-100-N00	Supporting Bar (Newton) 347-095-100-N00, Nominal value 5 N, stainless steel fine turned
347-115-100-N00	Supporting Bar (Newton) 347-115-100-N00, Nominal value 10 N, stainless steel fine turned

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-638N	DAkkS calibration certificate 962-638N
962-638RN	DAkkS Recalibration 962-638RN

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision are six days



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days