

Product brochure

08.08.2026

KERN 347-095-100

Supporting bar, OIML Class M1, 500 g, stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

FEATURES

- for fixing slotted weights
- The beam bar is tapered near the hook so that slotted weights can be hooked and unhooked easily. During use, the weights remain securely in place and are protected against accidentally falling off.
- Own weight beam bar 500 g
- Maximum total load 20 kg (is exclusive of the own weight of the beam bar, e.g. the maximum possible total weight is calculated from “Maximum total load” + “own weight beam bar”)
- Largest slotted weight possible 10 kg

PLEASE NOTE FOR YOUR ORDER

- Compatible with all slot weights/supporting bars of type A

PICTOGRAMS

Standard

Option





ACCESSORIES

Model	Description
347-015	Test weight 347-015 (OIML Class M1, Nominal value 1 g, Slot weight (Type A), stainless steel fine turned)
347-025	Test weight 347-025 (OIML Class M1, Nominal value 2 g, Slot weight (Type A), stainless steel fine turned)
347-035	Test weight 347-035 (OIML Class M1, Nominal value 5 g, Slot weight (Type A), stainless steel fine turned)
347-045	Test weight 347-045 (OIML Class M1, Nominal value 10 g, Slot weight (Type A), stainless steel fine turned)
347-055	Test weight 347-055 (OIML Class M1, Nominal value 20 g, Slot weight (Type A), stainless steel fine turned)
347-065	Test weight 347-065 (OIML Class M1, Nominal value 50 g, Slot weight (Type A), stainless steel fine turned)
347-075	Test weight 347-075 (OIML Class M1, Nominal value 100 g, Slot weight (Type A), stainless steel fine turned)
347-085	Test weight 347-085 (OIML Class M1, Nominal value 200 g, Slot weight (Type A), stainless steel fine turned)

Model	Description
347-095	Test weight 347-095 (OIML Class M1, Nominal value 500 g, Slot weight (Type A), stainless steel fine turned)
347-115	Test weight 347-115 (OIML Class M1, Nominal value 1 kg, Slot weight (Type A), stainless steel fine turned)
347-125	Test weight 347-125 (OIML Class M1, Nominal value 2 kg, Slot weight (Type A), stainless steel fine turned)
347-135	Test weight 347-135 (OIML Class M1, Nominal value 5 kg, Slot weight (Type A), stainless steel fine turned)
347-145	Test weight 347-145 (OIML Class M1, Nominal value 10 kg, Slot weight (Type A), stainless steel fine turned)
347-515-400	Plastic weight case 347-515-400 (for nominal values 500 g - 1 kg for classes M1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-639	DAkkS calibration certificate 962-639
962-639R	DAkkS Recalibration 962-639R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days