

Product brochure

08.08.2026

KERN 347-095-N00

Individual weight (Newton), 5 N, Slot weight (Type A), stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- Compatible with all slot weights/supporting bars of type A

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
317-280	Gloves 317-280 (cotton, 1 pair)
347-075-100-N00	Supporting Bar (Newton) 347-075-100-N00, Nominal value 1 N, stainless steel fine turned
347-095-100-N00	Supporting Bar (Newton) 347-095-100-N00, Nominal value 5 N, stainless steel fine turned
347-110-400	Plastic weight case 347-110-400 (for nominal values 1 kg for classes F1 - M3 for design Cylindrical/compact)
347-115-100-N00	Supporting Bar (Newton) 347-115-100-N00, Nominal value 10 N, stainless steel fine turned

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-639N	DAkkS calibration certificate 962-639N
962-639RN	DAkkS Recalibration 962-639RN

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision are six days



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days