

Product brochure

08.08.2026

KERN 347-1 15-100-N00

Supporting bar (Newton), 10 N, stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

FEATURES

- for fixing slotted weights
- The beam bar is tapered near the hook so that slotted weights can be hooked and unhooked easily. During use, the weights remain securely in place and are protected against accidentally falling off.

PLEASE NOTE FOR YOUR ORDER

- Compatible with all slot weights/supporting bars of type A

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
347-015-N00	Newton weight347-015-N00, Nominal value 1 cN, Slot weight (Type A), stainless steel fine turned
347-025-N00	Newton weight347-025-N00, Nominal value 2 cN, Slot weight (Type A), stainless steel fine turned
347-035-N00	Newton weight347-035-N00, Nominal value 5 cN, Slot weight (Type A), stainless steel fine turned
347-045-N00	Newton weight347-045-N00, Nominal value 10 cN, Slot weight (Type A), stainless steel fine turned
347-055-N00	Newton weight347-055-N00, Nominal value 20 cN, Slot weight (Type A), stainless steel fine turned
347-065-N00	Newton weight347-065-N00, Nominal value 50 cN, Slot weight (Type A), stainless steel fine turned
347-075-N00	Newton weight347-075-N00, Nominal value 1 N, Slot weight (Type A), stainless steel fine turned
347-085-N00	Newton weight347-085-N00, Nominal value 2 N, Slot weight (Type A), stainless steel fine turned
347-095-N00	Newton weight347-095-N00, Nominal value 5 N, Slot weight (Type A), stainless steel fine turned
347-115-N00	Newton weight347-115-N00, Nominal value 10 N, Slot weight (Type A), stainless steel fine turned
347-125-N00	Newton weight347-125-N00, Nominal value 20 N, Slot weight (Type A), stainless steel fine turned
347-135-N00	Newton weight347-135-N00, Nominal value 50 N, Slot weight (Type A), stainless steel fine turned
347-145-N00	Newton weight347-145-N00, Nominal value 100 N, Slot weight (Type A), stainless steel fine turned
347-515-400	Plastic weight case 347-515-400 (for nominal values 500 g - 1 kg for classes M1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-641N	DAkkS calibration certificate 962-641N
962-641RN	DAkkS Recalibration 962-641RN

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision are six days



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days