

Product brochure

08.08.2026

KERN 347-115

Individual weight, OIML Class M1, 1 kg, Slot weight (Type A), stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- Compatible with all slot weights/supporting bars of type A

PICTOGRAMS

Standard



Option



ACCESSORIES

Model	Description
317-280	Gloves 317-280 (cotton, 1 pair)
347-075-100	Supporting bar 347-075-100, OIML Class M1, Nominal value 100 g, stainless steel fine turned
347-095-100	Supporting bar 347-095-100, OIML Class M1, Nominal value 500 g, stainless steel fine turned
347-115-100	Supporting bar 347-115-100, OIML Class M1, Nominal value 1 kg, stainless steel fine turned
347-130-400	Plastic weight case 347-130-400 (for nominal values 5 kg for classes F1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-641	DAkks calibration certificate 962-641
962-641R	DAkks Recalibration 962-641R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days