

Product brochure

08.08.2026

# KERN 347-116

Individual weight, OIML Class M1, 1 kg, Hook weight, stainless steel fine turned

ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

+49 7433 9933-0  
info@kern-sohn.com  
www.kern-sohn.com

PICTOGRAMS

Standard

Option





ACCESSORIES

| Model       | Description                                                                                                      |
|-------------|------------------------------------------------------------------------------------------------------------------|
| 317-280     | Gloves 317-280 (cotton, 1 pair)                                                                                  |
| 347-120-400 | Plastic weight case 347-120-400 (for nominal values 2 kg   for classes F1 - M3   for design Cylindrical/compact) |

SERVICES

| Model    | Description                                                         |
|----------|---------------------------------------------------------------------|
| 920-350  | Marking by laser engraving (individual marking) KERN 920-350        |
| 920-351  | Marking of test weights (last digits of serial number) KERN 920-351 |
| 962-641  | DAkkS calibration certificate 962-641                               |
| 962-641R | DAkkS Recalibration 962-641R                                        |

ENGLISH

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days