

KERN 347-126

KERN

Individual weight, OIML Class M1, 2 kg, Hook weight, stainless steelfine turned



Shipping method	Parcel service
Net weight approx.	2,0 kg
Gross weight approx.	2,0 kg
Shipping weight	2 kg

Optional cases	
Case - Plastic article number	347-130-400

Pictograms

STANDARD



Category	
Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

Measuring System	
Test weight - volume (cm ³)	238,0952
Test weight - volume uncertainty (cm ³)	4,8186
Test weight - Density (kg/m ³)	8000
Test weight - Density uncertainty (kg/m ³)	100
Test weight - Density [Min] (OIML)	2600
Test weight - nominal value	2 kg
OIML-compliant	✓
Weight OIML class	M1
Permissible tolerance according to OIML (+/-)	100 mg

Construction	
Shape weight	Hook weight
Dimensions (Ø×H)	60×146 mm
Material	stainless steelfine turned

Services (optional)	
Article number for DAkkS calibration	962-642

Packing & Shipping	
Delivery time	1 d
Dimensions packaging (W×D×H)	60×60×146 mm