

Product brochure

08.08.2026

KERN 347-126-N00

Individual weight (Newton), 20 N, Hook weight, stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PICTOGRAMS

Standard

Option





ACCESSORIES

Model	Description
317-280	Gloves 317-280 (cotton, 1 pair)
347-130-400	Plastic weight case 347-130-400 (for nominal values 5 kg for classes F1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-642N	DAkkS calibration certificate 962-642N
962-642RN	DAkkS Recalibration 962-642RN

ENGLISH

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision are six days



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days