

Product brochure

08.08.2026

KERN 347-135-N00

Individual weight (Newton), 50 N, Slot weight (Type A), stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- Compatible with all slot weights/supporting bars of type A

PICTOGRAMS

Standard



Option



ACCESSORIES

Model	Description
317-280	Gloves 317-280 (cotton, 1 pair)
317-290	Gloves 317-290 (leather, cotton, 1 pair)
347-075-100-N00	Supporting Bar (Newton) 347-075-100-N00, Nominal value 1 N, stainless steel fine turned
347-095-100-N00	Supporting Bar (Newton) 347-095-100-N00, Nominal value 5 N, stainless steel fine turned
347-115-100-N00	Supporting Bar (Newton) 347-115-100-N00, Nominal value 10 N, stainless steel fine turned
347-140-400	Plastic weight case 347-140-400 (for nominal values 10 kg for classes F1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-643N	DAkKS calibration certificate 962-643N
962-643RN	DAkKS Recalibration 962-643RN

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision are six days



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days