

Product brochure

08.08.2026

# KERN 347-137-100

Supporting bar, OIML Class M1, 5 kg, stainless steel fine turned

ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

+49 7433 9933-0  
info@kern-sohn.com  
www.kern-sohn.com

## FEATURES

- for fixing slotted weights
- A precise curvature at the lower end of the beam bar ensures that the slotted weights rest securely. This design prevents the weights from slipping and ensures stability – even in slight inclined positions
- Own weight beam bar 5 kg

## PLEASE NOTE FOR YOUR ORDER

- Compatible with all slot weights/supporting bars of type B

## PICTOGRAMS

Standard

Option



## ACCESSORIES

| Model   | Description                                                                                                 |
|---------|-------------------------------------------------------------------------------------------------------------|
| 317-280 | Gloves 317-280 (cotton, 1 pair)                                                                             |
| 317-290 | Gloves 317-290 (leather, cotton, 1 pair)                                                                    |
| 347-037 | Test weight 347-037 (OIML Class M1, Nominal value 5 g, Slot weight (Type B), stainless steel fine turned)   |
| 347-047 | Test weight 347-047 (OIML Class M1, Nominal value 10 g, Slot weight (Type B), stainless steel fine turned)  |
| 347-057 | Test weight 347-057 (OIML Class M1, Nominal value 20 g, Slot weight (Type B), stainless steel fine turned)  |
| 347-067 | Test weight 347-067 (OIML Class M1, Nominal value 50 g, Slot weight (Type B), stainless steel fine turned)  |
| 347-077 | Test weight 347-077 (OIML Class M1, Nominal value 100 g, Slot weight (Type B), stainless steel fine turned) |
| 347-087 | Test weight 347-087 (OIML Class M1, Nominal value 200 g, Slot weight (Type B), stainless steel fine turned) |
| 347-097 | Test weight 347-097 (OIML Class M1, Nominal value 500 g, Slot weight (Type B), stainless steel fine turned) |
| 347-117 | Test weight 347-117 (OIML Class M1, Nominal value 1 kg, Slot weight (Type B), stainless steel fine turned)  |

| Model   | Description                                                                                                 |
|---------|-------------------------------------------------------------------------------------------------------------|
| 347-127 | Test weight 347-127 (OIML Class M1, Nominal value 2 kg, Slot weight (Type B), stainless steel fine turned)  |
| 347-137 | Test weight 347-137 (OIML Class M1, Nominal value 5 kg, Slot weight (Type B), stainless steel fine turned)  |
| 347-147 | Test weight 347-147 (OIML Class M1, Nominal value 10 kg, Slot weight (Type B), stainless steel fine turned) |

## SERVICES

| Model    | Description                                                         |
|----------|---------------------------------------------------------------------|
| 920-350  | Marking by laser engraving (individual marking) KERN 920-350        |
| 920-351  | Marking of test weights (last digits of serial number) KERN 920-351 |
| 962-643  | DAkkS calibration certificate 962-643                               |
| 962-643R | DAkkS Recalibration 962-643R                                        |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days