

Product brochure

08.08.2026

KERN 347-137

Individual weight, OIML Class M1, 5 kg, Slot weight (Type B), stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- Compatible with all slot weights/supporting bars of type B

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
317-280	Gloves 317-280 (cotton, 1 pair)
317-290	Gloves 317-290 (leather, cotton, 1 pair)
347-117-100	Supporting bar 347-117-100, OIML Class M1, Nominal value 1 kg, stainless steel fine turned
347-137-100	Supporting bar 347-137-100, OIML Class M1, Nominal value 5 kg, stainless steel fine turned
347-147-100	Supporting bar 347-147-100, OIML Class M1, Nominal value 10 kg, stainless steel fine turned

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-643	DAkks calibration certificate 962-643
962-643R	DAkks Recalibration 962-643R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days