

Product brochure

08.08.2026

KERN 347-146-N00

Individual weight (Newton), 100 N, Hook weight, stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

PICTOGRAMS

Standard

Option





ACCESSORIES

Model	Description
317-290	Gloves 317-290 (leather, cotton, 1 pair)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-644N	DAkkS calibration certificate 962-644N
962-644RN	DAkkS Recalibration 962-644RN

ENGLISH

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision are six days



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days