

Product brochure

08.08.2026

KERN 347-147-100

Supporting bar, OIML Class M1, 10 kg, stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

FEATURES

- for fixing slotted weights
- A precise curvature at the lower end of the beam bar ensures that the slotted weights rest securely. This design prevents the weights from slipping and ensures stability – even in slight inclined positions
- Own weight beam bar 10 kg

PLEASE NOTE FOR YOUR ORDER

- Compatible with all slot weights/supporting bars of type B

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
317-290	Gloves 317-290 (leather, cotton, 1 pair)
347-037	Test weight 347-037 (OIML Class M1, Nominal value 5 g, Slot weight (Type B), stainless steel fine turned)
347-047	Test weight 347-047 (OIML Class M1, Nominal value 10 g, Slot weight (Type B), stainless steel fine turned)
347-057	Test weight 347-057 (OIML Class M1, Nominal value 20 g, Slot weight (Type B), stainless steel fine turned)
347-067	Test weight 347-067 (OIML Class M1, Nominal value 50 g, Slot weight (Type B), stainless steel fine turned)
347-077	Test weight 347-077 (OIML Class M1, Nominal value 100 g, Slot weight (Type B), stainless steel fine turned)
347-087	Test weight 347-087 (OIML Class M1, Nominal value 200 g, Slot weight (Type B), stainless steel fine turned)
347-097	Test weight 347-097 (OIML Class M1, Nominal value 500 g, Slot weight (Type B), stainless steel fine turned)
347-117	Test weight 347-117 (OIML Class M1, Nominal value 1 kg, Slot weight (Type B), stainless steel fine turned)
347-127	Test weight 347-127 (OIML Class M1, Nominal value 2 kg, Slot weight (Type B), stainless steel fine turned)

Model	Description
347-137	Test weight 347-137 (OIML Class M1, Nominal value 5 kg, Slot weight (Type B), stainless steel fine turned)
347-147	Test weight 347-147 (OIML Class M1, Nominal value 10 kg, Slot weight (Type B), stainless steel fine turned)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-644	DAkKS calibration certificate 962-644
962-644R	DAkKS Recalibration 962-644R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days