

Product brochure

08.08.2026

KERN 347-147-N00

Individual weight (Newton), 100 N, Slot weight (Type B), stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- Compatible with all slot weights/supporting bars of type B

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
317-290	Gloves 317-290 (leather, cotton, 1 pair)
347-117-100-N00	Supporting Bar (Newton) 347-117-100-N00, Nominal value 10 N, stainless steel fine turned
347-137-100-N00	Supporting Bar (Newton) 347-137-100-N00, Nominal value 50 N, stainless steel fine turned
347-147-100-N00	Supporting Bar (Newton) 347-147-100-N00, Nominal value 100 N, stainless steel fine turned

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-644N	DAkks calibration certificate 962-644N
962-644RN	DAkks Recalibration 962-644RN

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision are six days



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days