


Product brochure

08.08.2026

KERN 348-03

Milligram weights, M1, 5 mg, Platelet, stainless steel polished

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
315-247	Tweezers 315-247 (stainless steel, plastic 136 mm)
317-009-600	Aluminium weight case 317-009-600 (for nominal values 1 mg - 500 mg (einzel) for classes E1 - M3 for design Fracture gram)
318-270	Brush 318-270
338-090-200	Wood weight case 338-090-200 (for nominal values 1 mg - 500 mg (einzel) for classes E1 - M2 for design Cylindrical/compact)
347-009-400	Plastic weight case 347-009-400 (for nominal values 1 mg - 500 mg (einzel) for classes E1 - M2 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
962-653	DAkks calibration certificate 962-653
962-653R	DAkks Recalibration 962-653R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days