


KERN 348-09

KERN

Milligram weights, M1, 500 mg, Platelet, stainless steel polished



Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Optional cases

Case - Plastic article number	347-009-400
Case - Aluminium article number	317-009-600
Case - Wood article number	338-090-200

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	348

Measuring System

Conformity assessment possible	✓
Test weight - Density (kg/m ³)	7960
Test weight - Density uncertainty (kg/m ³)	170
Test weight - nominal value	500 mg
OIML-compliant	✓
Weight OIML class	M1
Permissible tolerance according to OIML (+/-)	0,8 mg

Construction

Shape weight	Platelet
Dimensions weight (W×D×H)	12×12×4 mm
Material	stainless steel polished

Services (optional)

Article number for conformity assessment	952-659
Article number for DAkkS calibration	962-659

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	12×12×4 mm
Shipping method	Parcel service