


Product brochure

08.08.2026

# KERN 356-89

Individual weight, OIML Class M2, 50 kg, Block, cast iron lacquered

ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

+49 7433 9933-0  
info@kern-sohn.com  
www.kern-sohn.com

PICTOGRAMS

Standard

Option





ACCESSORIES

| Model   | Description                              |
|---------|------------------------------------------|
| 317-290 | Gloves 317-290 (leather, cotton, 1 pair) |

SERVICES

| Model    | Description                                                  |
|----------|--------------------------------------------------------------|
| 920-350  | Marking by laser engraving (individual marking) KERN 920-350 |
| 962-646  | DAkkS calibration certificate 962-646                        |
| 962-646R | DAkkS Recalibration 962-646R                                 |

ENGLISH

## KERN Pictogram Overview



Pallet shipment: The duration of internal product provision are two days



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days