

# KERN 357-01

Individual weight, OIML Class M2, 1 g, Knob, stainless steelfine turned



Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

## Optional cases

Case - Plastic article number	347-030-400
Case - Aluminium article number	317-010-600
Case - Wood article number	337-010-200

## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML M2
Product family	357

## Measuring System

Test weight - volume (cm <sup>3</sup> )	0,119
Test weight - volume uncertainty (cm <sup>3</sup> )	0,0024
Test weight - Density (kg/m <sup>3</sup> )	7950
Test weight - Density uncertainty (kg/m <sup>3</sup> )	140
Test weight - nominal value	1 g
OIML-compliant	✓
Weight OIML class	M2
Permissible tolerance according to OIML (+/-)	3 mg

## Construction

Shape weight	Knob
Dimensions (Ø×H)	6×6 mm
Material	stainless steelfine turned

## Services (optional)

Article number for DAkkS calibration	962-631
--------------------------------------	---------

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	24×24×51 mm
Shipping method	Parcel service