

KERN 357-11

KERN

Individual weight, OIML Class M2, 1 kg, Knob, stainless steelfine turned



Shipping method	Parcel service
Net weight approx.	1,00 kg
Gross weight approx.	1,00 kg
Shipping weight	1 kg

Optional cases	
Case - Plastic article number	347-110-400
Case - Aluminium article number	317-110-600
Case - Wood article number	337-110-200

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M2
Product family	357

Measuring System

Test weight - volume (cm ³)	119,0476
Test weight - volume uncertainty (cm ³)	2,4093
Test weight - Density (kg/m ³)	7950
Test weight - Density uncertainty (kg/m ³)	140
Test weight - Density [Min] (OIML)	2300
Test weight - nominal value	1 kg
OIML-compliant	✓
Weight OIML class	M2
Permissible tolerance according to OIML (+/-)	160 mg

Construction

Shape weight	Knob
Dimensions (Ø×H)	48×80 mm
Material	stainless steelfine turned

Services (optional)

Article number for DAkkS calibration	962-641
--------------------------------------	---------

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	61×61×94 mm