

Product brochure

08.08.2026

KERN 366-92

Individual weight, OIML Class M3, 200 g, Cylinders, cast iron lacquered

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

PICTOGRAMS

Standard

Option





SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
962-638	DAkkS calibration certificate 962-638
962-638R	DAkkS Recalibration 962-638R

ENGLISH

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days