

Product brochure

08.08.2026

KERN 367-07

Individual weight, OIML Class M3, 100 g, Knob, stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
315-250	Tweezers 315-250 (stainless steel 150 mm)
317-070-600	Aluminium weight case 317-070-600 (for nominal values 100 g for classes E1 - M3 for design Cylindrical/compact)
337-070-200	Wood weight case 337-070-200 (for nominal values 100 g for classes F2 + M1 for design Cylindrical/compact)
347-070-400	Plastic weight case 347-070-400 (for nominal values 50 g - 100 g for classes F1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-637	DAkkS calibration certificate 962-637
962-637R	DAkkS Recalibration 962-637R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days