

Product brochure

08.08.2026

KERN 367-09

Individual weight, OIML Class M3, 500 g, Knob, stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
315-245	Tweezers 315-245 (stainless steel 250 mm)
317-090-600	Aluminium weight case 317-090-600 (for nominal values 500 g for classes E1 - M3 for design Cylindrical/compact)
317-280	Gloves 317-280 (cotton, 1 pair)
337-090-200	Wood weight case 337-090-200 (for nominal values 500 g for classes F2 + M1 for design Cylindrical/compact)
347-090-400	Plastic weight case 347-090-400 (for nominal values 500 g for classes F1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-639	DAkks calibration certificate 962-639
962-639R	DAkks Recalibration 962-639R

ENGLISH

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days