

Product brochure

08.08.2026

# KERN 367-12

Individual weight, OIML Class M3, 2 kg, Knob, stainless steel fine turned

ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

+49 7433 9933-0  
info@kern-sohn.com  
www.kern-sohn.com

## PICTOGRAMS

Standard

Option



## ACCESSORIES

| Model       | Description                                                                                                        |
|-------------|--------------------------------------------------------------------------------------------------------------------|
| 315-245     | Tweezers 315-245 (stainless steel   250 mm)                                                                        |
| 317-120-600 | Aluminium weight case 317-120-600 (for nominal values 2 kg   for classes E1 - M3   for design Cylindrical/compact) |
| 317-280     | Gloves 317-280 (cotton, 1 pair)                                                                                    |
| 337-120-200 | Wood weight case 337-120-200 (for nominal values 2 kg   for classes F2 + M1   for design Cylindrical/compact)      |
| 347-120-400 | Plastic weight case 347-120-400 (for nominal values 2 kg   for classes F1 - M3   for design Cylindrical/compact)   |

## SERVICES

| Model    | Description                                                         |
|----------|---------------------------------------------------------------------|
| 920-350  | Marking by laser engraving (individual marking) KERN 920-350        |
| 920-351  | Marking of test weights (last digits of serial number) KERN 920-351 |
| 962-642  | DAkks calibration certificate 962-642                               |
| 962-642R | DAkks Recalibration 962-642R                                        |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days