

KERN OZB-A45 15UK-2026e

KERN

Double gooseneck illumination for best illumination of your sample, with plug-in power supply for UK



Gross weight approx.	2,0 kg
Shipping weight	3,8 kg
Shipping volume weight	3,8 kg

Product information

GTIN/EAN number	404576 1490826
Trade number	94054239
REACH classification	no prohibited substances

Pictograms

STANDARD



Category

Brand	Optics
Product category	Microscope
Product group	Microscope illumination
Product family	OZB-IF

Construction

Dimension (W×D×H)	120×370×250 mm
Gooseneck length	300 mm
Beleuchtungsart	LED
Illuminance	6 W
Colour temperature	5600-6300 K
Illumination dimmable	✓

Power Supply

Supplied power supply	Power supply unit
Plug-in power supply / adapter for countries - optional	UK

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Approval

CE mark	✓
---------	---

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	300×250×250 mm
Shipping method	Parcel service
Net weight approx.	1,6 kg