

Equipment qualification			
Logbook for qualification of balance and weighing systems			
General details of equipment			
Device Type:	Single-range balance <input checked="" type="checkbox"/>	Multi-range balance <input type="checkbox"/>	Multi-interval balance <input type="checkbox"/>
Verified Weighing Device:	<input type="checkbox"/>		
	Range 1	Range 2	Range 3
Weighting range (Max):	320 g	N/A	N/A
Scale interval (d):	0,001 g	N/A	N/A
Verification value (d):	0,001 g	N/A	N/A
Min:	0,01 g	N/A	N/A
Connected accessories:	Printer YKG-01		
Operating instructions managed by:	Process center		
Date of installation:	28-08-2023		
Date of putting out of service:	N/A		
Service contract number:	N/A		
Service case telephone number:	N/A		
Equipment qualification executed by:	<div style="border: 1px solid black; padding: 2px;"> qualification officer </div>		
Time			
Name of qualification officer:	<div style="border: 1px solid black; padding: 2px;"> (signature) (date) </div>		
Equipment qualification reviewed by:	<div style="border: 1px solid black; padding: 2px;"> (signature) (date) </div>		
John Doe	<div style="border: 1px solid black; padding: 2px;"> (signature) </div>		
Name of measuring and test equipment officer:			

Erstellungsdatum: 23.06.2026

Kategorie

Marke	KERN
Produktkategorie	KERN Software
Produktgruppe	Geräte-Qualifizierung
Produktfamilie	SEQ

Display

Sprachen des Benutzerinterfaces	Deutsch Englisch Französisch Spanisch Italienisch
---------------------------------	---

Verpackung & Versand

Lieferzeit	1 d
Versandart	Download

Produktinformationen

GTIN/EAN-Nummer	4045761361720
Warennummer (Zolltarifnummer)	85234910
REACH-Einstufung	keine Stoffe über der zulässigen Konzentration

Piktogramme