

Equipment qualification			
Logbook for qualification of balance and weighing systems			
General details of equipment			
Device Type:	<input checked="" type="checkbox"/> Single-range balance	<input type="checkbox"/> Multi-range balance	<input type="checkbox"/> Multi-interval balance
Verified Weighing Device:	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Weighting range (Max):	Range 1 320 g	Range 2 N/A	Range 3 N/A
Scale interval (d):	0.001 g	N/A	N/A
Verification value (g):	0.001 g	N/A	N/A
Min:	0.01 g	N/A	N/A
Connected accessories:	Printer YKG-01		
Operating instructions managed by:	Process owner		
Date of installation:	28-08-2023		
Date of putting out of service:	N/A		
Service contract number:	N/A		
Service case telephone number:	N/A		
Equipment qualification executed by:	<div style="border: 1px solid black; padding: 2px;"> <div style="text-align: right; font-size: small;">qualification officer</div>  </div>		
Time & other:			
(Name of qualification officer):	<div style="border: 1px solid black; padding: 2px;"> <div style="text-align: right; font-size: small;">(signature)</div>  </div>		
Equipment qualification reviewed by:	<div style="border: 1px solid black; padding: 2px;"> <div style="text-align: right; font-size: small;">EQS</div> </div>		
John Doe	<div style="border: 1px solid black; padding: 2px;"> <div style="text-align: right; font-size: small;">(signature)</div>  </div>		
(Name of measuring and test equipment officer):			

Data di creazione: 23.06.2026

Categoria

Marchio	KERN
Categoria di prodotto	KERN Software
Gruppo di prodotti	Qualificazione dello strumento
Famiglia di prodotti	SEQ

Display

Linguaggi dell'interfaccia utente	Tedesco Inglese Francese Spagnolo Italiano
-----------------------------------	--

Imballaggio e spedizione

Tempi di consegna	1 d
Metodo di spedizione	Download

Informazioni sul prodotto

GTIN/Numero EAN	4045761361720
Numero merceologico	85234910
Classificazione REACH	nessuna sostanza proibita

Pittogrammi