

SAUTER TVP-L

KERN

Manual test stand for compressive force measurements with digital length measurement



Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	570×275×240 mm
Shipping method	Parcel service
Net weight approx.	11 kg
Gross weight approx.	12 kg
Shipping weight	12 kg

Pictograms

STANDARD



Category

Brand	Sauter
Product category	Test bench
Product group	Manual test stand
Product family	TVP-L

Measuring System

Measurable length [Max]	10 cm
Tolerance (measured length)	2%

Construction

Dimension (W×D×H)	150×233×420 mm
Dimension housing (W×D×H)	150×233×420 mm
Dimensions completely mounted (W×D×H)	150×233×420 mm
Material housing	metal, lacquered
Type of drive	manual lever
Mounting - force dissipation	5 × threaded hole M6
Mounting - force application	5 × threaded hole M6
Generated force [Max]	500 N
Measuring room height [Max]	315 mm
Traverse distance with lever movement / release [Max]	78 mm

Environmental conditions

Storage temperature [Min]	-10 °C
Storage temperature [Max]	40 °C

Approval

CE mark	✓
---------	---