

KERN YDB-02

KERN

Set for density determination of liquids and solids with density $\leq/\geq 1$



Pictograms

STANDARD



Category

Brand	KERN
Product category	Balance
Product group	Kit for density determination
Product family	YDB

Construction

Material	stainless steel
----------	-----------------

Environmental conditions

Storage temperature [Min]	-10 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	300×445×185 mm
Shipping method	Parcel service
Net weight approx.	1,6 kg
Gross weight approx.	2,6 kg
Shipping weight	4,9 kg
Shipping volume weight	4,9 kg

Product information

GTIN/EAN number	4045761504455
Trade number	90160090
REACH classification	no prohibited substances