

# KERN YDB-03

Universal density set

**KERN**



## Pictograms

STANDARD



## Category

Brand	KERN
Product category	Balance
Product group	Kit for density determination
Product family	YDB

## Construction

Material	stainless steel
----------	-----------------

## Environmental conditions

Storage temperature [Min]	-10 °C
Storage temperature [Max]	40 °C

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	390×245×210 mm
Shipping method	Parcel service
Net weight approx.	0,65 kg
Gross weight approx.	1,6 kg
Shipping weight	4 kg
Shipping volume weight	4 kg

## Product information

GTIN/EAN number	4045761504462
Trade number	90160090
REACH classification	no prohibited substances